



**UNIVERSITY OF COLORADO  
HOSPITAL**

ANSCHUTZ MEDICAL CAMPUS

**Cytomegalovirus IGG (CMVIGG) and  
Cytomegalovirus IGM (CMVIGM)**

Effective on April 15, 2009 the Clinical Laboratory will begin Cytomegalovirus IGM and Cytomegalovirus IGG testing by Chemiluminescent Immunoassay.

	<b>Cytomegalovirus Antibody IgG</b>	<b>Cytomegalovirus Antibody IgM</b>
TEST CODE	CMVIGG	CMVIGM
TEST METHOD	Chemiluminescent Immunoassay	Chemiluminescent Immunoassay
PERFORMED	TUESDAY	TUESDAY
REPORTED	Same day	Same day
SPECIMEN	Serum	Serum
MIN VOLUME	0.5 mL	0.5 mL
STABILITY	Serum refrigerated, 2 days Serum frozen (-20C), 1 month	Serum refrigerated, 2 days Serum frozen (-20C), 1 month
REFERENCE RANGE	Negative	Negative
CPT CODE	86644	86645

Please call Mary Ellen Boyle, Chief Technologist at 720-848-7067 if you have any questions or visit our website at <http://www.uch.edu/for-healthcare-professional/Clinical-Laboratory/index.aspx> for additional information.

Issued [4-16-09]

**Clinical Laboratory**  
**Medical Director Ronald B. Lepoff, MD**