



**University of Colorado Hospital**

ANSCHUTZ MEDICAL CAMPUS

## **Test Update AMIKACIN Method Change**

The immunoassay test method for serum Amikacin on the TDx analyzer has been discontinued by the manufacturer. Starting on Wednesday, November 30, 2011, serum Amikacin testing will be performed by an immunoturbidimetric testing method on the Beckman DXC analyzer. The new test will be available on demand 24 hours per day. Reference ranges will remain the same. There will be no change to current CPT code, sample type, or order names.

Please call Mary Ellen Boyle at 720-848-7067 if you have any questions or visit our website at <http://www.uch.edu/for-healthcare-professional/Clinical-Laboratory/index.aspx> for additional information.

Issued November 28, 2011

**Clinical Laboratory**  
**Medical Director Ronald B. Lepoff, MD**