



**UNIVERSITY OF COLORADO**  
**HOSPITAL**

ANSCHUTZ MEDICAL CAMPUS

**Test Update**  
**Vitamin D (VITD25)**

On February 2, 2009 University of Colorado Hospital, Clinical Laboratory will begin testing for the quantitative determination of 25-hydroxy Vitamin D and other hydroxylated Vitamin D metabolites in serum.

<b>Test mnemonic</b>	<b>VITD25</b>
Method	Chemiluminescent Immunoassay
Sample	Serum
Volume	0.500 mL
Minimum Volume	0.300 mL
Storage/Stability	Serum, refrigerated (2-8 °C), 5 days, Serum, frozen (-20 °C), 1 month
CPT	82306
Performed	M, W, F
Reported	Same day

Please call Lorne Holland, MD at 720-848-7050 if you have any questions or visit our website at <http://www.uch.edu/for-healthcare-professional/Clinical-Laboratory/index.aspx> for additional information.

Issued [2-2-2009]

**Clinical Laboratory**  
**Medical Director Ronald B. Lepoff, MD**