



UNIVERSITY OF COLORADO
HOSPITAL

ANSCHUTZ MEDICAL CAMPUS

Test Update
Partial Thromboplastin Time (PTT) Reference Range

The University Of Colorado Hospital Clinical Laboratory has revised the reference ranges for the following tests to reflect the sensitivity of the test methods:

TEST	OLD Reference Range	NEW Reference Range
PTT	23.8 – 33.5 seconds	24.9 – 36.0 seconds
Therapeutic Range for Heparin Monitoring	60.7 – 88.1 seconds	62.5 – 91.1 seconds

Please call Stuart Lind, MD at 303-724-1383 if you have any questions or visit our website at <http://www.uch.edu/for-healthcare-professional/Clinical-Laboratory/index.aspx> for additional information.

Issued [11-03-2009]

Clinical Laboratory
Medical Director Ronald B. Lepoff, MD