

Weekly Liquid QC Log									
Date	Initials	Level #	Lot #	Exp. Date	THb Result	THb Acceptable Range	O2Hb Result	O2Hb Acceptable Range	Cuvette Pathlength
		Level 1							
		Level 2							
		Level 3							

Daily Optical QC Log	
Yellow Filter Serial Num: _____	Orange Filter Serial Number: _____
Yellow THb Range: _____	Orange THb Range: _____
Yellow O2Hb Range: _____	Orange O2Hb Range: _____

Date	Time	Initials	Yellow Filter		Orange Filter	
			O2Hb Result	THb Result	O2Hb Result	THb Result
Day 1						
Day 2						
Day 3						
Day 4						
Day 5						
Day 6						
Day 7						

*Liquid QC is performed with each new lot number, new shipment of cuvettes, and once a week. Optical QC is performed every 8 hours of patient testing.

User Initials	Print Name	User Initials	Print Name	User Initials	Print Name